

# STAINLESS STEEL LABORATORY CARTS, BASKETS AND TRAYS

MODEL LC-1004 / LC-1003D / LTP1 / LTP2 / STP1 / STP2



Model LC-1004



Model LC-1003D

## Applications

The Tuff stainless steel laboratory carts are designed to be used in high workload teaching, training and research laboratories.

The carts are tailored to support the material handling needs of the popular dryers that are being marketed by our company throughout the years. A series of stainless steel baskets, trays and test tube racks are designed to accommodate drying stocks, samples or glasswares of different sizes and configurations for use in the laboratories.

The baskets and trays fit well on the laboratory carts and provide laboratory technicians the ease of transporting the drying stocks, glasswares, samples or apparatus from the laboratory to the cleaning area, drying room, or equipment store or from laboratory to laboratory.

## Stainless Steel Reversible Laboratory Carts

Tuff offers two models of laboratory carts. The model LC-1004 has four shelves while the Model LC-1003D has three shelves and two equipment drawers. Both stainless steel carts can withstand chemical spills and can be steam-cleaned.

These carts can be moved safely and easily. Both carts have 19 gauge stainless steel shelves and legs are made of 1" OD, 18 gauge, 304 series hollow stainless steel tubing complete with 3" heavy duty swivel castors (two with press lock). Both the models have 3" edges on the sides and back of all shelves. Convenient handles on both sides of carts reduce handle strain by keeping hands at comfortable level.

## Specifications

| MODEL / CAT. NO. | DIMENSIONS              | NO. OF SHELVES | DRAWERS | SHELF CLEARANCE | LOAD   | WEIGHT |
|------------------|-------------------------|----------------|---------|-----------------|--------|--------|
| LC-1004          | 1040 W x 600 D x 1040 H | 4 units        | -       | 250 mm          | 200 kg | 37 kg  |
| LC-1003D         | 1020 W x 580 D x 1040 H | 3 units        | 2 units | 300 mm          | 300 kg | 46 kg  |

| STANDARD PACKAGE FOR MODEL FDD-720 / FDD-1000 / FD-825 |   |         |            |
|--|---|---------|------------|
| Model  | Description   | Picture | Cat. No.   |
| LTP1   | 4-shelf Laboratory Cart complete with 8 units of empty Type A Basket  |         | LC-1004-8A |
| LTP2   | 4-shelf Laboratory Cart complete with 8 units of empty Type H Tray of dimensions 527 W x 377 D x 20 H mm each |         | LC-1004-8H |

| STANDARD PACKAGE FOR MODEL FSD-380 / FD-380 DRYERS |  |         |             |
|--|--|---------|-------------|
| Model  | Description  | Picture | Cat. No.    |
| STP1   | 3-shelf 2-drawer Laboratory Cart complete with 6 units of empty Type A Basket  |         | LC-1003D-6A |
| STP2   | 3-shelf 2-drawer Laboratory Cart complete with 6 units of empty Type H Tray of dimensions 527 W x 377 D x 20 H mm each |         | LC-1003D-6H |

All specifications are subject to change without prior notice.