



Model LC-1004



Model LC-1003D

Applications

The Tuff stainless steel laboratory carts are designed to be used in high workload teaching, training and research laboratories.

The carts are tailored to support the material handling needs of the popular dryers that are being marketed by our company throughout the years. A series of stainless steel baskets, trays and test tube racks are designed to accommodate drying stocks, samples or glasswares of different sizes and configurations for use in the laboratories.

The baskets and trays fit well on the laboratory carts and provide laboratory technicians the ease of transporting the drying stocks, glasswares, samples or apparatus from the laboratory to the cleaning area, drying room, or equipment store or from laboratory to laboratory.

Stainless Steel Reversible Laboratory Carts

Tuff offers two models of laboratory carts. The model LC-1004 has four shelves while the Model LC-1003D has three shelves and two equipment drawers. Both stainless steel carts can withstand chemical spills and can be steam-cleaned.

These carts can be moved safely and easily. Both carts have 19 gauge stainless steel shelves and legs are made of 1" OD, 18 gauge, 304 series hollow stainless steel tubing complete with 3" heavy duty swivel castors (two with press lock). Both the models have 3" edges on the sides and back of all shelves. Convenient handles on both sides of carts reduce handle strain by keeping hands at comfortable level.

Specifications

MODEL / CAT. NO.	DIMENSIONS (mm)	NO. OF SHELVES	DRAWERS	SHELF CLEARANCE	LOAD	WEIGHT
LC-1004	1040 W x 600 D x 1040 H	4 units	-	250 mm	200 kg	37 kg
LC-1003D	1020 W x 580 D x 1040 H	3 units	2 units	300 mm	300 kg	46 kg

STANDARD PACKAGE FOR MODEL FDD-720 / FDD-1000 / FD-825			
Model	Description	Picture	Cat. No.
LTP1	4-shelf Laboratory Cart complete with 8 units of empty Type A Basket		LC-1004-8A
LTP2	4-shelf Laboratory Cart complete with 8 units of empty Type H Tray of dimensions 527 W x 377 D x 20 H mm each		LC-1004-8H

STANDARD PACKAGE FOR MODEL FSD-380 / FD-380 DRYERS			
Model	Description	Picture	Cat. No.
STP1	3-shelf 2-drawer Laboratory Cart complete with 6 units of empty Type A Basket		LC-1003D-6A
STP2	3-shelf 2-drawer Laboratory Cart complete with 6 units of empty Type H Tray of dimensions 527 W x 377 D x 20 H mm each		LC-1003D-6H

All specifications are subject to change without prior notice.