



Single-Door Cabinets Specifications

MODEL NO.	EXTERNAL DIMENSIONS W x D x H mm(in.)	INTERNAL DIMENSIONS W x D x H mm(in.)	CAPACITY Gal. (L)	No. of Shelves	Door Type
CCG-004-AA-A	435 x 435 x 630 (17 x 17 x 25)	352 x 352 x 481 (14 x 14 x 19)	4(15)	1	Self-close
CCG-012-AA-A	595 x 460 x 960 (23 x 18 x 38)	512 x 377 x 811 (20 x 15 x 32)	12(45)	1	Self-close
CCG-015-AA-A	595 x 460 x 1190 (23 x 18 x 47)	512 x 377 x 1041 (20 x 15 x 41)	15(57)	1	Self-close
CCG-004-BA-A	435 x 435 x 630 (17 x 17 x 25)	352 x 352 x 481 (14 x 14 x 19)	4(15)	1	Manual
CCG-012-BA-A	595 x 460 x 960 (23 x 18 x 38)	512 x 377 x 811 (20 x 15 x 32)	12(45)	1	Manual
CCG-015-BA-A	595 x 460 x 1190 (23 x 18 x 47)	512 x 377 x 1041 (20 x 15 x 41)	15(57)	1	Manual
CCH-022-AA-A	595 x 460 x 1720 (23 x 18 x 68)	512 x 378 x 1574 (20 x 15 x 62)	22(83)	3	Self-close
CCH-054-AA-A	595 x 860 x 1720 (23 x 34 x 68)	512 x 778 x 1574 (20 x 31 x 62)	54(204)	3	Self-close
CCH-022-BA-A	595 x 460 x 1720 (23 x 18 x 68)	512 x 378 x 1574 (20 x 15 x 62)	22(83)	3	Manual
CCH-054-BA-A	595 x 860 x 1720 (23 x 34 x 68)	512 x 778 x 1574 (20 x 31 x 62)	54(204)	3	Manual



Double-Door Cabinets Specifications

MODEL NO.	EXTERNAL DIMENSIONS W x D x H mm(in.)	INTERNAL DIMENSIONS W x D x H mm(in.)	CAPACITY Gal. (L)	No. of Shelves	Door Type
CCE-030-AA-A	1090 x 460 x 1190 (43 x 18 x 47)	1007 x 375 x 1041 (40 x 15 x 41)	30(114)	1	Self-close
CCE-045-AA-A	1090 x 460 x 1720 (43 x 18 x 68)	1007 x 375 x 1574 (40 x 15 x 62)	45(170)	2	Self-close
CCE-060-AA-A	860 x 860 x 1720 (34 x 34 x 68)	778 x 778 x 1574 (31 x 31 x 62)	60(227)	2	Self-close
CCE-090-AA-A	1090 x 860 x 1720 (43 x 34 x 68)	1007 x 778 x 1574 (40 x 31 x 62)	90(340)	2	Self-close
CCE-030-BA-A	1090 x 460 x 1190 (43 x 18 x 47)	1007 x 375 x 1041 (40 x 15 x 41)	30(114)	1	Manual
CCE-045-BA-A	1090 x 460 x 1720 (43 x 18 x 68)	1007 x 375 x 1574 (40 x 15 x 62)	45(170)	2	Manual
CCE-060-BA-A	860 x 860 x 1720 (34 x 34 x 68)	778 x 778 x 1574 (31 x 31 x 62)	60(227)	2	Manual
CCE-090-BA-A	1090 x 860 x 1720 (43 x 34 x 68)	1007 x 778 x 1574 (40 x 31 x 62)	90(340)	2	Manual

All specifications are subject to change without prior notice.