

FLAMMABLE SAFETY STORAGE CABINETS
SLIDING-DOOR CABINETS
 MODEL CCF SERIES



Sliding-Door Cabinets
Specifications



MODEL NO.	EXTERNAL DIMENSIONS W x D x H mm(in.)	INTERNAL DIMENSIONS W x D x H mm(in.)	CAPACITY Gal. (L)	No. of Shelves	Door Type
CCF-030-AA-A	1090 x 460 x 1190 (43 x 18 x 47)	1007 x 375 x 1041 (40 x 15 x 41)	30(114)	1	Right hand / Self-close
CCF-045-AA-A	1090 x 460 x 1720 (43 x 18 x 68)	1007 x 375 x 1574 (40 x 15 x 62)	45(170)	2	Right hand / Self-close
CCF-060-AA-A	860 x 860 x 1720 (34 x 34 x 68)	778 x 778 x 1574 (31 x 31 x 62)	60(227)	2	Right hand / Self-close
CCF-090-AA-A	1090 x 860 x 1720 (43 x 34 x 68)	1007 x 778 x 1574 (40 x 31 x 62)	90(340)	2	Right hand / Self-close
CCF-030-BA-A	1090 x 460 x 1190 (43 x 18 x 47)	1007 x 375 x 1041 (40 x 15 x 41)	30(114)	1	Right hand / Manual
CCF-045-BA-A	1090 x 460 x 1720 (43 x 18 x 68)	1007 x 375 x 1574 (40 x 15 x 62)	45(170)	2	Right hand / Manual
CCF-060-BA-A	860 x 860 x 1720 (34 x 34 x 68)	778 x 778 x 1574 (31 x 31 x 62)	60(227)	2	Right hand / Manual
CCF-090-BA-A	1090 x 860 x 1720 (43 x 34 x 68)	1007 x 778 x 1574 (40 x 31 x 62)	90(340)	2	Right hand / Manual

All specifications are subject to change without prior notice.

LOW CORROSIVE LIQUID STORAGE CABINETS
SINGLE-DOOR CABINETS & DOUBLE-DOOR CABINETS
 MODEL CCG / CCE SERIES



Specifications

MODEL NO.	EXTERNAL DIMENSIONS W x D x H mm(in.)	INTERNAL DIMENSIONS W x D x H mm(in.)	CAPACITY Gal. (L)	No. of Shelves	Door Type
CCG-012-DA-B	595 x 460 x 960 (23 x 18 x 38)	512 x 377 x 811 (20 x 15 x 32)	12(45)	1	Single door / Manual
CCG-015-DA-B	595 x 460 x 1190 (23 x 18 x 47)	512 x 377 x 1041 (20 x 15 x 41)	15(57)	1	Single door / Manual
CCE-030-BA-B	1090 x 460 x 1190 (43 x 18 x 47)	1007 x 375 x 1041 (40 x 15 x 41)	30(114)	1	Double door / Manual
CCE-045-BA-B	1090 x 460 x 1720 (43 x 18 x 68)	1007 x 375 x 1574 (40 x 15 x 62)	45(170)	2	Double door / Manual
CCE-060-BA-B	860 x 860 x 1720 (34 x 34 x 68)	778 x 778 x 1574 (31 x 31 x 62)	60(227)	2	Double door / Manual
CCE-090-BA-B	1090 x 860 x 1720 (43 x 34 x 68)	1007 x 778 x 1574 (40 x 31 x 62)	90(340)	2	Double door / Manual

All specifications are subject to change without prior notice.

FLAMMABLE SAFETY STORAGE CABINETS / LOW CORROSIVE LIQUID STORAGE CABINETS