



Model CSC-9500



Model CSC-91000



Model CSC-91000C

### Features

- Total elimination of noxious vapours by ventilation fan.
- Low running cost as no expensive consumables are required e.g. sampling kits, test kits or chemical filters.
- Safe to use as there is no risk of non-functioning saturated filters.
- No hassle as there is no need to have routine or scheduled samplings or checks on filters.
- Constructed of heavy gauge mild steel plate coated with anti-corrosion powder paint.
- Doors constructed of PMMA which does not support burning in case of fire.
- Excellent visual comfort with large doors which allows full view of the stored chemicals.
- Key lock security prevents access to unauthorised persons.
- Magnetic labelling tags.
- “Fan on” indicator light.
- Adjustable fan speed.
- Silent operation.
- Adjustable shelves – Four supplied for each storage compartment.
- Supplied with spillage tray(s) and 3 meters of flexible duct piping.

### Specifications

MODEL	CSC-9500	CSC-9500C	CSC-91000	CSC-91000C
Type	Single compartment Single door	Double compartment Double-door	Double compartment Double-door	Four compartment Four-door
Capacity	100 one-litre bottles and flasks	100 one-litre bottles and flasks	200 one-litre bottles and flasks	200 one-litre bottles and flasks
Applications	For the storage of acids, bases or solvents	For the storage of acids, bases or acid-solvents or solvents-acids in 2 isolated compartments	For the storage of acids, bases or solvents in 2 isolated compartments	For the storage of acids, bases or solvents in 4 isolated compartments
Ventilation rate	Up to 275 m <sup>3</sup> /hr	Up to 275 m <sup>3</sup> /hr	Up to 550 m <sup>3</sup> /hr	Up to 550 m <sup>3</sup> /hr
Fan speed	Adjustable	Adjustable	Adjustable	Adjustable
Outlet duct diameter	4"	4"	4"	4"
Shelf dimensions (mm)	697W x 420D x 20H (4 units)	320W x 420D x 20H (8 units)	718W x 415D x 20H (8 units)	330W x 415D x 20H (16 units)
Internal dimensions of each compartment (mm)	720W x 470D x 1690H	345W x 470D x 1690H (Two compartments)	745W x 470D x 1710H (Two compartments)	355W x 470D x 1710H (Four compartments)
External dimensions (mm)	800W x 605D x 2125H	800W x 605D x 2125H	1600W x 605D x 2150H	1600W x 605D x 2150H
Spillage tray capacity	20 litres	20 litres	20 litres (2 units)	20 litres (2 units)
Electrical supply	240 VAC, 50 Hz, 0.4 A	240 VAC, 50 Hz, 0.4 A	240 VAC, 50 Hz, 0.8 A	240 VAC, 50 Hz, 0.8 A
Weight	145 kg	155 kg	325 kg	345 kg

All specifications are subject to change without prior notice.