



Chemstore Bench-Type / Wall-Mount / Under-Counter Flammable Liquid Storage Cabinets For Limited Work / Floor Space

These highly visible and durable yellow powder-coated cabinets feature a “HIGHLY FLAMMABLE – KEEP FIRE AWAY” warning on the front doors. These storage cabinets are constructed of double-walled, 18-gauge steel with 1½-inch (3.8 cm) insulating air barrier. Features include piano hinges, three-point lock with handle, dual vents with flash arrestors, 2-inch (5 cm) raised, leakproof door sill to contain spills and four levelling legs. Four brackets support each adjustable shelf. Each shelf can support up to 350 lbs (160 kg) load.

Specifications

MODEL	CS-320	CS-317	CS-322
Type	Bench-type	Wall-mount	Under-counter
Capacity	20 gallons	17 gallons	22 gallons
No. of shelves supplied	1 unit	1 unit	1 unit
Shelf dimensions (mm)	992W x 215D x 35H	992W x 370D x 35H	805W x 470D x 35H
Interior dimensions (mm)	1012W x 240D x 1040H	1012W x 380D x 480H	820W x 480D x 790H
Exterior dimensions (mm)	1090W x 305D x 1140H	1092W x 460D x 610H	900W x 560D x 900H
Weight	80.6 kg	61 kg	74.3 kg



Model CS-317



Model CS-322

Chemstore Bench-Type / Wall-Mount / Under-Counter Acid and Corrosive Liquid Storage Cabinets For Limited Work / Floor Space

- 18-gauge steel cabinet
- Double-wall construction
- Superior powder-coated paint finish
- Dual vents
- Raised leak-proof door sill
- Three-point latching doors with lock and handle
- Full-length piano hinges
- Adjustable levelling feet
- Reinforced fibre glass spill-catcher tray shelves supplied as standard
- Tray shelves can be easily removed and cleaned
- Reinforced adjustable tray shelves to support up to 350 lbs (160 kg) load each
- Not recommended for storage of phenol, nitric acid or sulfuric acid.

Specifications

MODEL	ASC-320	ASC-317	ASC-322
Type	Bench-type	Wall-mount	Under-counter
Capacity	250 litres	180 litres	310 litres
No. of tray shelves supplied	2 units	2 units	2 units
Shelf dimensions (mm)	992W x 215D x 35H	992W x 370D x 35H	805W x 470D x 35H
Interior dimensions (mm)	1012W x 240D x 1040H	1012W x 380D x 480H	820W x 480D x 790H
Exterior dimensions (mm)	1090W x 305D x 1140H	1092W x 460D x 610H	900W x 560D x 900H
Weight	83.6 kg	63 kg	77.7 kg

All specifications are subject to change without prior notice.

MANUFACTURED IN ACCORDANCE WITH OSHA REGULATIONS AND NFPA CODE 30 SPECIFICATIONS