



Pharmaceutical Refrigerators

Up to 3 or 8 units of refrigerators



Blood Bank Refrigerators



Hi-Low Scientific Refrigerators



GSM Remote Alarm System



### Stand Alone GSM Remote Alarm System

- Monitors up to three or eight units of medical refrigerators simultaneously.
- Dials up to four mobile numbers in case of alarm conditions.

#### Dial-up conditions:

- When any medical refrigerator alarm is triggered.
- When GSM control module is in battery-low condition.
- When GSM system is armed or unarmed.

#### Note:

- The system works independently from the main temperature controller of the medical refrigerators.
- All Mediref Medical Refrigerators and Protech Hi-low Scientific Refrigerators (old or new units) can be upgraded to the new system.
- Customers can buy the system with an option of one, two or three channels, depending on their requirements / budgets



CO<sub>2</sub> Back-up System



Chart Recorder

### CO<sub>2</sub> Back-up System

The CO<sub>2</sub> Back-up System is supplied with a battery back-up system and all hardware except the liquid CO<sub>2</sub> cylinder and gas regulator.

### Chart Recorder

All Mediref Medical Refrigerators if needed, can be supplied with a 7-day chart recorder with a temperature span of -35°C to 15°C. The scale is 1mm/°C and the accuracy is ±2°C.