



TSC-400G



TSC-900G



TSC-900GS



TSC-900GM

To cater for laboratory storage requirements, we now manufacture a wide range of specially designed metal storage cabinets to suit your needs.

These cabinets come in modular form as free standing, wall-hung or mobile units.

Free Standing Metal Storage Cabinets:

The free standing cabinets are available in single or double door, full glass or half glass, solid door or open shelves units. You may have choice of hinged door or sliding door models. All hinged doors are provided with silent pads to reduce noise during closure.

The cabinets and shelves are constructed of 1.2mm mild steel plate powder coated with epoxy paint for durability. To meet different applications, you may have the choice of stainless steel shelves, extra stainless steel tray or fiber glass tray where hygiene, spillage or resistant to chemicals or solvent is a consideration for your application. The shelves are full length adjustable to accommodate storage of different product types.

Specifications

Type	Model	Description	Width (mm)	Depth (mm)	Height (mm)
Free Standing Metal Storage Cabinets	TSC-400M	Solid door	500	480	1800
	TSC-400GM	Half glass door	500	480	1800
	TSC-400G	Glass door	500	480	1800
	TSC-900M	Solid door	1145	480	1800
	TSC-900GM	Half glass door	1145	480	1800
	TSC-900G	Glass door	1145	480	1800
	TSC-900GS	Glass sliding door	1145	480	1800
	TSC-900	Open shelves	1145	445	1800

All specifications are subject to change without prior notice.

Optional Accessories

Cabinet	Standard Shelf	Stainless steel shelf	Stainless steel tray	Fiberglass tray
Model 400	C/N TSC-400S	C/N TSC-400SS	C/N TSC-400T	C/N TSC-400FT
Model 900	C/N TSC-900S	C/N TSC-900SS	C/N TSC-900T	C/N TSC-900FT

All specifications are subject to change without prior notice.