



Model S-230SS



Model S-460SS

The solution warming cabinets are available in three models. One model is bench-top cabinet with capacity of 230 litres (Model S-230SS). The other model is a floor standing model with 230 litres warming cabinet at the top compartment and 230 litres of storage space at the bottom compartment (Model S-460SS). Also we are offering another model with two warming compartments (model S-460SS-2C)

Applications

The Solution Warming Cabinet is used to warm and maintain temperature of infusion solutions, irrigation fluids, filled syringes, ointment lotions, etc.

Features

- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set temperature at all times
- Built-in timer
- Uniform heating
- Fully insulated on all side walls
- Audible and visual over-temperature alarm
- Glass or stainless steel door(s)



Specifications

Model	S-230SS	S-460SS	S-460SS-2C
Type	Bench-top, single-door	Floor-standing, double-door	
Working Capacity	200 litres	200 litres and 200 litres storage	400 litres
Preset temperature	42°C	42°C (top compartment only)	42°C
Working temperature range	35°C - 50°C	35°C - 50°C (top compartment only)	35°C - 50°C
Safety devices	Built-in thermostat and RCCB. Fitted with temperature alarm to avoid overheating.		
Timer	Built-in 999 minutes timer		
No. of shelves supplied	2 units	4 units (2x2)	4 units (2x2)
Shelf material	Stainless steel		
Shelf dimensions (mm)	485W x 425D x 15H		
Interior dimensions (mm)	540W x 480D x 800H	540W x 480D x 800H (Two chambers of 200 litres each)	540W x 480D x 800H (Two chambers of 200 litres each)
Exterior dimensions (mm)	800W x 640D x 1030H	800W x 640D x 1915H	800W x 640D x 1915H
Power consumption	1000 watt	1000 watt	2000 watt
Electrical	220-240 VAC, 50Hz, 5A	220-240 VAC, 50Hz, 5A	220-240 VAC, 50Hz, 9A
Weight	66 kg	114 kg	120 kg

All specifications are subject to change without prior notice.