



Model HP-3002D



Model HP-3003D

**Features**

- Extra large hot plate designed for service laboratories and educational laboratories.
- Large work surface to accommodate a large number of test containers / samples.

**Applications**

- For preheating of a large volume of test samples.
- For general heating applications, dissolving of chemicals or media, etc.
- For digestion of a large volume of samples.

**Specifications**

| MODEL                     | HP-3002D                           | HP-3003D              | HP-3004D           |
|---------------------------|------------------------------------|-----------------------|--------------------|
| Size of hot plate (feet)  | 2.0W x 1.5D                        | 3.0W x 1.5D           | 4.0W x 2.0D        |
| Power                     | 2.0 kW                             | 3.0 kW                | 4.0 kW             |
| Temperature               | Up to 300°C                        | Up to 200°C           | Up to 200°C        |
| Temperature control       | Digital controller, display to 1°C |                       |                    |
| Electrical                | 220-240 VAC, 50Hz                  | 415 VAC, 50Hz, 3φ, 4W |                    |
| Overall Dimensions (feet) | 2.0W x 1.5D x 0.5H                 | 3.0W x 1.5D x 0.5H    | 4.0W x 2.0D x 0.5H |
| Weight                    | 18 kg                              | 27 kg                 | 35 kg              |

**HOTPLATE WITH OPTIONAL SAND BATH ATTACHMENT**  
 MODEL HP-3001DS



Model HP-3001DS

**Features**

- Constructed of Stainless steel.
- Bath height (side panel) 2".
- With mounting base that sits/secure the bath on the hot plate.

**Specifications**

| MODEL                     | HP-3001DS               |
|---------------------------|-------------------------|
| Size of hot plate (mm)    | 430 W x 290 D x 25 H    |
| Power                     | 2.0 kW                  |
| Temperature               | Ambient to 300°C        |
| Temperature control       | Digital controller      |
| Electrical                | 220-240 VAC, 50Hz, 7.5A |
| Overall Dimensions (feet) | 435 W x 340 D x 155 H   |
| Weight                    | 12 kg                   |

All specifications are subject to change without prior notice.