



**Features**

The Tuff stainless steel laboratory carts and baskets will enhance the productivity and workflow in the high workload laboratories in teaching, training and research institutes and production laboratories.

You can efficiently load and transport the complete drying stock in one go from the cleaning area to the drying room and to the laboratory and vice versa.



The tray baskets can be used as draining racks or as holding or transport trays in the laboratory.



The tray baskets hold the glasswares vertically effectively increasing the capacity of the dryer.

**Specifications of Laboratory Carts**

MODEL/CAT.	DIMENSION	SHELVES	DRAWERS	SHELF CLEARANCE	LOAD
LC-1004	1040 W x 600 D x 1040 H	4 units	-	250 mm	200 kg
LC-1003D	1020 W x 580 D x 1040 H	3 units	2 units	300 mm	300 kg

**Ordering Information (Standard Packages)**

For Model FDD-720 / FDD-1000

For Model FSD-380



Model : LTP1  
C/N: LC-1004-8A  
Laboratory Cart with 8 units Type A Basket



Model : LTP2  
C/N: LC-1004-8H  
Laboratory Cart with 8 units Type H Tray



Model : STP1  
C/N: LC-1003D-6A  
Laboratory Cart with 6 units Type A Basket



Model : STP2  
C/N: LC-1003D-6H  
Laboratory Cart with 6 units Type H Tray

**Ordering Information (Standard Items)**

ITEM	DESCRIPTION	CAT. NO.
Model LC-1004 Laboratory Cart	Laboratory cart with 4 shelves (can hold 8 baskets or trays)	LC-1004
Model LC-1003D Laboratory Cart	Laboratory cart with 3 shelves and 2 drawers (can hold 6 baskets or trays)	LC-1003D
Type A Basket	Empty basket for loading samples (Dimensions: 460 W x 478 D x 170 H mm)	SSB
Type B Basket	Basket with perforated tray (121 holes of 30 x 30 mm each)	SSB-3030
Type C Basket	Basket with perforated tray (42 holes of 50 x 50 mm each)	SSB-5050
Type D Basket	Basket with perforated tray (21 holes of 50 x 50 mm & 60 holes of 30 x 30 mm)	SSB-3050
Type E Basket	Basket with 6 units of Z-Shaped Type S test tube rack (40 holes of 13 mm Ø each)	SSB-6S
Type F Basket	Basket with 6 units of Z-Shaped Type M test tube rack (40 holes of 16 mm Ø each)	SSB-6M
Type G Basket	Basket with 6 units of Z-Shaped Type L test tube rack (24 holes of 20 mm Ø each)	SSB-6L
Type H Tray	Stainless steel tray of dimensions 527 W x 377 D x 20 H mm	SST
Z-Shaped Test Tube Rack (Type S)	Test tube rack with 40 holes of 13 mm Ø each	TTRM-4013
Z-Shaped Test Tube Rack (Type M)	Test tube rack with 40 holes of 16 mm Ø each	TTRM-4016
Z-Shaped Test Tube Rack (Type L)	Test tube rack with 24 holes of 20 mm Ø each	TTRM-2420

All specifications are subject to change without prior notice.