

# STAINLESS STEEL LABORATORY CARTS and BASKETS

Designed for use with Protech High Performance Dryers



The Tuff stainless steel laboratory carts have the capacity to hold eight stainless steel baskets or trays (for Model LC-1004) or six stainless steel baskets or trays (for Model LC-1003D).

The baskets are constructed of perforated stainless steel plates that allow efficient ventilation and do not obstruct hot air flow from reaching the samples.



## Specifications of Laboratory Carts

MODEL/CAT.	DIMENSION	SHELVES	DRAWERS	SHELF CLEARANCE	LOAD
LC-1004	1040 W x 600 D x 1040 H	4 units	-	250 mm	200 kg
LC-1003D	1020 W x 580 D x 1040 H	3 units	2 units	300 mm	300 kg

## Ordering Information (Standard Packages)

### For Model FD-825



Model : LTP1  
C/N: LC-1004-8A  
Laboratory Cart with 8 units Type A Basket



Model : LTP2  
C/N: LC-1004-8H  
Laboratory Cart with 8 units Type H Tray

### For Model FD-380



Model : STP1  
C/N: LC-1003D-6A  
Laboratory Cart with 6 units Type A Basket



Model : STP2  
C/N: LC-1003D-6H  
Laboratory Cart with 6 units Type H Tray

## Ordering Information (Standard Items)

ITEM / DESCRIPTION	MODEL / CATALOGUE NO.
Laboratory Cart with 4 shelves	LC-1004
Laboratory Cart with 3 shelves and 2 drawers	LC-1003D
Type A Basket	SSB
Type H Tray	SST

All specifications are subject to change without prior notice.