



Model DW-FL90 / DW-FL270



Temperature controller

Features

- Vertical cabinet (upright)
- Operating temperature ranging from -10°C to -40°C can be set freely
- Microprocessor temperature controller
- LCD display
- Door with lock
- Audible / visual alarm : High or low temperature alarm, system failure alarm
- Danfoss compressor and EBM fan motor
- 4 casters provided for easy mobility
- Refrigerant : R507, CFC free
- Certificates : CE mark for DW-FL90 / DW-FL270, ISO9001, ISO14001, ISO13485
- OPTIONAL : 7- day chart recorder

Specifications

Model	DW-FL90	DW-FL270
CE MARK	Yes	Yes
Capacity	90 litres	270 litres
Drawers/shelves (fixed position)	4 drawers	7 drawers
Exterior construction	Powder coated mild – steel exterior	
Cabinet construction	ABS interior	
Insulation	90 mm polyurethane foam	
External size (mm)	623W x 544D x 842H	623W x 544D x 1680H
Internal size (mm)	438W x 388D x 539H	438W x 388D x 1580H
Input power	280 watt	370 watt
Electrical	220 – 240 VAC, 50 Hz, single phase	
Weight	44/48 kg	65/70 kg

All specifications are subject to change without prior notice.