



Features and Specifications

- Upright type
- Inside temperature ranging from -40°C to -86°C can be set freely
- Stainless steel inner door
- Supplied with access port as standard
- Heated door seal prevents ice build-up
- Efficient and well distributed evaporator with maximum contact with walls provides fast cool down
- Stainless steel interior, exterior is powder coated steel panel
- Front opening lockable door with full height handle
- 120mm foaming insulation with VIP insulation board
- Vacuum-release port
- Inner 3 shelves, two compartments with stainless steel inner door (DW-HL218 and DW-HL100 are 2 shelves, three compartments)
- Digital temperature display
- Keyboard lock and password protected configuration page
- Audible and visual alarm system:
High or low temperature alarm, Power failure alarm, Low battery alarm, Door open alarm, Filter blocking alarm, System failure alarm
- 72 hour battery backup for printer and controller display
- 4 units casters for easy mobility
- Certificates: ISO9001, ISO13485
- Power supply: 220-240 VAC, 50 Hz, 1 phase
- Optional: CO₂ back up, storage racks, boxes, temperature monitoring system

Specifications

Model	Capacity	External size (W x D x H) mm	Internal size (W x D x H) mm	Input Power	No. of inner door	No. of Shelf	Weight (Net / Gross)
DW-HL100 (CE Mark)	100 litres	720 x 793 x 1000	450 x 450 x 500	800W	2	1	120/140KG
DW-HL218 (CE Mark)	218 litres	940 x 840 x 1552	582 x 470 x 768	1000W	2	1	225/245KG
DW-HL398S (CE Mark)	398 litres	785 x 1041 x 1947	440 x 696 x 1266	1000W	2	3	300/320KG
DW-HL528S (CE Mark)	528 litres	930 x 1041 x 1947	585 x 696 x 1266	1100W	2	3	320/340KG
DW-HL678S (CE Mark)	678 litres	1090 x 998 x 1968	750 x 696 x 1268	1200W	4	3	340/360KG
DW-HL778S (CE Mark)	778 litres	1205 x 998 x 1968	865 x 696 x 1268	1200W	4	3	350/370KG
DW-HL858S (CE Mark)	858 litres	1218 x 997 x 1998	877 x 696 x 1378	1300W	2	2	440/460KG
DW-HL1008S (CE Mark)	1008 litres	1363 x 997 x 1998	1022 x 696 x 1378	1300W	2	2	480/500KG